

# Liquid Bacterial Culture

## Overview

Grow bacteria such as *E. coli* or *B. subtilis* in nutrient-rich broth, allowing rapid multiplication in a suspended, aerated environment.

### A. Materials

- a. Plate with bacterial colonies
- b. Liquid culture tubes
- c. TB medium
- d. Antibiotics
  - i. Chloramphenicol ( )
  - ii. Ampicillin (25mg/ml)
  - iii. Erythromycin ( )
- e. Bunsen burner
- f. Pipette
- g. Incubator

### B. Procedure

1. Label liquid culture tubes
2. Add 3 mL of TB medium and antibiotics at a 1:1,000 dilution to each culture tube.
3. Dip a pipette on the selected colony
4. Dispense the tip into tube
5. Place the tube in the incubator and incubate for 16 hours with the tube positioned at a slant angle.

## Liquid Culture

